

①

M. HAUSCHILD
CERN
26-Oct-00

dE/dx Particle ID at the TESLA-TPC

- estimated dE/dx resolution
(based on running detectors experience)
- toy Monte Carlo results
(study by Magali Gruwe)
- particle separation power

②

Do we need Particle ID?

- Physics studies show

- only "weak" physics cases, e.g.
background reduction at heavy flavor tags
 - no need for dedicated particle ID detector
- dE/dx of TPC sufficient ("for free")
not compromising other detector performances

- Different situation at "Giga Z" option

= $10^9 Z^0$ at $2 \times 45 \text{ GeV}$

- Z^0 physics largely profits from particle ID
e.g. CP violation studies in b decays
(Richard Hawkings)



reconstruction of $B^0 \rightarrow \pi^+ \pi^-$ decays

no dE/dx

$S/B = 10:1$

with dE/dx (moderate res.)

$S/B = 40:1$



3)

Estimated dE/dx Resolution

(based on running detectors experience)

- Lehraus 1983:


$$\frac{\sigma dE/dx}{dE/dx} = 5.7\% * L^{-0.37}$$

↗ "total detector depth"
= track length * pressure

- with present detectors:

$$... = 5.5\% * L^{-0.36} \quad \leftarrow \text{very similar}$$

- do interpolation between DELPHI ($= 5.7\%$)
ALEPH ($= 4.5\%$)




④


Estimated dE/dx Resolution II

cccccccccccccccccccc

- Slightly better correlation using
 - total number of electrons (primary + secondary)
 - number of clusters (primary electrons only)



← TESLA - TPC :
 Ar/CH_4 (90:10)
 $120 \text{ cm track length}$
 1 bar




5

Toy Monte Carlo Studies

ooooooooooooooo

- More detailed studies needed (toy MC)
 - Various theoretical models on dE/dx available:
 - e.g. Allison, Cobb
 - Sternheimer, Peierls
 - Bethe, Landau ...
 - not always consistent picture
+ some complicated + difficult to use
- try to avoid any dE/dx model
 - use measured cluster size distributions
(Fischle, Heintze, Schmidt NIM A301, 102, 1031)
 - generate # of clusters along track (Poisson)
+ generate cluster size according
to measured distributions



Measured Cluster Size Distributions

(also available: He, CO₂)

⑥


Toy Monte Carlo

- Add:

- diffusion ($\sim 1\text{m}$)
 - gas gain fluctuations (10^4)
 - electronics noise
 - crosstalk between samples } Vary + check influence
- Correlations between samples ↗

- Get Landau distributions

Compare with measured ones (here: OPAL)




7

Landau Distributions

ooooooooooooooo

- Keep track length constant ($120\text{ cm} = \text{TESLA-TPC}$)
 - Vary sampling length (or number of samples)

Landau distributions, 600 e⁻ noise, Cross-talk

↖ ALEPH pre-amps

- More (shorter) samples → better resolution ?

$$\text{e.g. } \sigma(dE/dx) \sim N^{-0.43} ?$$

↖ Not true in general ! ↗ valid for constant sampling length only

↖ Many empty samples if sampling length too short


↖ difficult to handle in later dE/dx reconstruction

→ Sampling length needs optimization

(8)


Truncated Mean

- Use standard truncated mean method \nwarrow get rid of to reconstruct dE/dx per track (max. likelihood might be better but less robust) \nwarrow nowhere used at any running detector
- Resolution depends on upper truncation (where to cut the Landau tail?)



optimal truncation
depends on
sampling length


best resolution
at 75% - 80%
truncation for
100 ... 300 samples



9

dE/dx Resolution (opt.)


- Saturation reached at $\sim 200 - 300$ samples
($\sim 4 - 5$ mm sampling length)



- Again: Do not use too many samples (too short sampling length)

reconstructed dE/dx
becomes very sensitive
to hit detection
threshold for
 > 300 samples
(< 4 mm sampling length)



BIAS!



10

Simulation \leftrightarrow Measurements

	simulated	measured
ALEPH	4.2 %	4.5 %
DELPHI	5.2 %	5.7 %
TOPAZ	3.5 %	4.4 %
TPC/2 γ	2.5 %	3.0 %
OPAL	2.6 %	2.8 %
TESLA	4.3 % with 240 samples (0.5 cm) 4.6 % with 120 samples (1 cm)	



⑪

Particle Separation Power

- Important for physics:

$$\text{particle separation power} = \frac{(dE/dx)_A - (dE/dx)_B}{5(dE/dx)_{A,B}}$$

in relativistic rise ! ($> 2 \text{ GeV}$)

✓ separation
↑ resolution

- dE/dx separation depends on


- gas mixture:

	ionization (resolution)	rel. rise (separation)
light gas (He etc.)	bad	large ✓
hydrocarbons	good ✓	small

- gas pressure:

high pressure

better resolution
→ smaller rel. rise






better resolution
dominates, best
separation power
at 3-4 bar

← 1.58

↑
but high pressure
detector not
feasible

⑫ Expected Separation Power

- Simulation + measured dE/dx resolutions based on
 - isolated, clean tracks, full number of samples
 - "test beam like"
- Particle ID in physics analysis use
 - tracks in dense track environments (multihadronic jets)
 - reduced no. of samples (limited by double hit res.)
 - samples less clean (additional corrections)
 - at ALEPH, DELPHI, OPAL typically
 $\sim 60 - 70\%$ of full no. of samples used



dE/dx resolutions
(incl. 0.2-0.3 % calibration error)

- - - 4.8 % (4.6 % w/o cal. error)
(TESLA-CDR, 120 samples)
- - - 4.3 % (4.1 % w/o cal. error)
(TESLA, 240 samples + opt. trunc.)
- - - 6.1 % (5.8 % w/o cal. error)
(TESLA, dense track env.)

4.3 % best
6.1 % average

dE/dx particle separation power

(13)

Conclusions

Conclusion

- dE/dx in TPC only tool
for (hadronic) particle ID at TESLA
- results of detailed toy MC study (Magali Gruwe)
 - able to reproduce measured detector resolutions
(with small offsets, calibration ?)
 - optimal resolution with $\sim 5\text{ mm}$ sampling length
(240 samples)
 - some improvements possible using optimized truncation
- estimated separation power
(based on 6.1% average dE/dx resolution
in multihadrons)
 - $e/\pi > 25$ up to 15.2 GeV
 - $\pi/K > 25$ up to 13.6 GeV
 - $K/p \sim 1.25$ max. in rel. rise
- Continue studies
 - different gases
 - influence of pad structure
 - (mixture of different sampling lengths)
at GEM pads
 - double hits / dense track environment
 - realistic dE/dx reconstruction